

Datasheet

ZNF37A monoclonal antibody (M01), clone 5F9

Catalog Number: H00007587-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant ZNF37A.

Clone Name: 5F9

Immunogen: ZNF37A (NP_001007095, 121 a.a. ~ 230 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

NKNGNSFWLNEDLIWHQKIKNWEQSFEYNECGKAFP
ENSLFLVHKRGYTGQKTCKYTEHKGKTCDSFFITHQQ
THPRENHYGNECGENIFEESILLEHQSVPFSQKLNL

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.2

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 7587

Gene Symbol: ZNF37A

Gene Alias: FLJ3472, KOX21, ZNF37