

## Datasheet

### KIF13B monoclonal antibody (M01), clone 6E11

**Catalog Number:** H00023303-M01

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against a partial recombinant KIF13B.

**Clone Name:** 6E11

**Immunogen:** KIF13B (NP\_056069, 1 a.a. ~ 100 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

**Sequence:**

MGDSKVKVAVRIRPMNRRETLHTKCVVDVDANKVIL  
NPVNTNLSKGDARGQPKVFAYDHCFSMDESVEKEY  
AGQDIVFKCLGENILQNAFDGYNACIF

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, S-ELISA, WB-Ce, WB-Re  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Isotype:** IgG2a Kappa

**Storage Buffer:** In 1x PBS, pH 7.2

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 23303

**Gene Symbol:** KIF13B

**Gene Alias:** GAKIN, KIAA0639